


Drinking Water Semivolatiles at 10ng on Rxi™-1ms

- | | | |
|--|-----------------------------------|---------------------------------------|
| 1. 2-fluorophenol (surr.) | 21. 2-naphthalenamine | 42. metolachlor |
| 2. bis(2-chloroethyl)ether | 22. 5-nitro- <i>o</i> -toluidine | 43. fluoranthene |
| 3. phenol-d6 (surr.) | 23. diethylphthalate | 44. pyrene |
| 4. 1,4-dichlorobenzene-d4 (int. std.)• | 24. fluorene | 45. butachlor |
| 5. nitrobenzene-d5 (surr.) | 25. propachlor | 46. <i>p</i> -terphenyl-d14 (surr.) |
| 6. naphthalene-d8 (int. std.)• | 26. diphenylamine | 47. <i>p</i> -dimethylaminoazobenzene |
| 7. naphthalene | 27. 2,4,6-tribromophenol (surr.) | 48. benzyl butyl phthalate |
| 8. 1-methylnaphthalene | 28. simazine | 49. 2-acetylaminofluorene |
| 9. 2-methylnaphthalene | 29. prometon | 50. bis(2-ethylhexyl)adipate |
| 10. hexachlorocyclopentadiene | 30. atrazine | 51. benzo(a)anthracene |
| 11. EPTC | 31. hexachlorobenzene | 52. chrysene-d12 (int. std.)• |
| 12. 2-fluorobiphenyl (surr.) | 32. 4-aminobiphenyl | 53. chrysene |
| 13. 2,6-dinitrotoluene | 33. terbacil | 54. bis(2-ethylhexyl)phthalate |
| 14. dimethylphthalate | 34. phenanthrene-d10 (int. std.)• | 55. di- <i>n</i> -octylphthalate |
| 15. acenaphthylene | 35. phenanthrene | 56. benzo(b)fluoranthene |
| 16. acenaphthene-d10 (int. std.)• | 36. anthracene | 57. benzo(k)fluoranthene |
| 17. acenaphthene | 37. metribuzin | 58. benzo(a)pyrene |
| 18. 2,4-dinitrotoluene | 38. acetochlor | 59. perylene-d12 (int. std.)• |
| 19. 1-naphthalenamine | 39. alachlor | 60. indeno(1,2,3-cd)pyrene |
| 20. molinate | 40. bromacil | 61. dibenzo(a,h)anthracene |
| | 41. di- <i>n</i> -butylphthalate | 62. benzo(ghi)perylene |



Column: Rxi™-1ms, 30m, 0.25mm ID, 0.25µm (cat.# 13323)
Sample: US EPA Method 525.2 mix: custom 525.2 calibration mix, SV Internal Standard Mix (cat.# 31206), B/N Surrogate Mix (4/89 SOW) (cat.# 31024), Acid Surrogate Mix (4/89 SOW) (cat.# 31025)
Inj.: 1.0µL, 10µg/mL each analyte (internal standards 100µg/mL), split (10:1) 4mm Drilled Uniliner® inlet liner (hole at bottom) (cat.# 20756)
Instrument: Agilent 6890
Inj. temp.: 250°C
Carrier gas: helium, constant flow
Flow rate: 1.2mL/min.
Oven temp.: 50°C (hold 1 min.) to 265°C @ 20°C/min., to 330°C @ 6°C/min. (hold 1 min.)
Det.: Agilent 5973 MSD
Transfer line temp.: 280°C
Scan range: 35-550 amu
Solvent delay: 3.20 min.
Tune: DFTPP
Ionization: EI

• Internal standards at 100ng on-column.